# Study of Mode Characteristics of Glass Optical Fibers LABORATORY MANUAL: TESTER MCS-04

- 1. Definition and Computation of Numerical Aperture
- 2. Definition and Computation of V Number & Computation
- 3. Specifications of 9/125 Single Mode SI Glass Fiber
- 4. Specifications of 50/125 Multimode GI Glass Fiber
- 5. Specifications of 62.5/125 Multimode GI Glass Fiber
- 6. Contents of Tester MCS04
- 7. Full Set-up to Study Single Mode Fiber Characteristics
- 8. Close-up of Laser Diode Unit
- 9. Close-up of Detector Screen
- 10. LP01 Intensity Pattern Observed on the Screen
- 11. LP11 Intensity Pattern Observed on the Screen
- 12. LP11 Intensity Pattern Observed on the Screen
- 13. LP22 Intensity Pattern Observed on the Screen
- 14. NA Computation from Screen Shot of LP01 Pattern
- 15. Full Set-up to Study Multimode GI Fiber Characteristics
- 16. Close-up of Laser Diode Unit
- 17. Close-up of Detector Screen Diode Unit
- 18. Intensity Patterns Observed on the Screen
- 19. NA Computation from Screen Shot of Multimode Fiber Spot

TELENET SYSTEMS, 203, Mehtab Arcade, TARNAKA, Secunderabad 500017 Telangana (info@telenet-systems.com) **1. Definition and Computation of Numerical Aperture** 

Numerical Aperure of an optical fiber is the measure of its Acceptance Angle.

Acceptance Angle is defined as the maximum angle at which the optical fiber accepts and transmits light through its core.

Numerical Aperture (NA) is also defined by the equation given below:

where  $n_o = 1$  (refractive index of air),  $\theta_c$  is half the Acceptance Angle of the fiber.

Numerical Aperture (NA) is also defined by the refractive indices of the optical fiber as given below:

 $NA = (n_{core}^2 - n_{cladding}^2)^{1/2}$ 

where  $n_{core}$  and  $n_{cladding}$  are the refractive indices of the fiber core and the fiber cladding

Qualitatively, NA is a measure of the light gathering ability of a fiber.

A number of optical parameters of the optical fiber are dependent on NA.

### 2. Definition and Computation of V Number & N

The normalized frequency, V of a fiber (also called the V number or the cut-off frequency) is given by the equation:

V=2∏a/ (NA)

where a is the core radius,  $\lambda$  is the wavelength of light, NA is the numerical aperture and  $n_{core}$  &  $n_{cladding}$  are the refractive indices of the fiber core and fiber cladding respectively

V number is also related to the number of guided modes, N, in a fiber.

For a step index fiber N =  $V^2/2$ .

For a graded index fiber (with a parabolic refractive index core profile),  $N = V^2/4$ .

If V is less than 2.405 for a given wavelength, the fiber is single mode.

If V is greater than 2.405 for a given wavelength, the fiber is multimode.

3. Specifications of 9/125 Single Mode SI Glass Fiber

## **Fiber Type**

**Numerical Aperture Core Diameter** Cladding Diameter 125+/- 1 microns Secondary Coating 250 microns Attenuation @ 1310 nm 0.34 db per Km Attenuation @ 1550 nm 0.20 db per Km **Effective Group Refractive Index** Cut- off wavelength Mode Field Diameter @ 9.2 +/- 0.4 microns 1310 nm Mode Field Diameter @ 10.4 +/- 0.6 microns 1550 nm **Refractive Index** Difference

**Glass Single Mode Step Index (9/125)** 0.11 8.3 microns 1.4675 @ 1310 nm 1.4681 @ 1550 nm 1260 nm 0.36%

4. Specifications of 50/125 Multimode GI Glass Fiber

# **Fiber Type**

Operating Wavelength Numerical Aperture **Core Diameter** Cladding Diameter 125+/- 1 microns Secondary Coating 250 microns Effective Group Refractive Index

Glass Multi Mode Graded Index (50/125) 750-1450 nm 0.20+/- 0.015 50 +/- 2.5 microns Bandwidth @ 850nm Greater than 500 MHz-Km Bandwidth @ 1300 nm Greater than 1000 MHz-Km Attenuation @ 850 nm Better than 2.8 db per Km Attenuation @ 1300 nm Better than 0.9 db per Km 1.480 @ 850nm 1.479 @ 1300nm

5. Specifications of 62.5/125 Multimode GI Glass Fiber

# Fiber Type

Operating Wavelength Numerical Aperture **Core Diameter** Cladding Diameter Secondary Coating 250 microns Attenuation @ 850 nm Attenuation @ 1300 nm **Effective Group** Refractive Index


**Glass Multi Mode** Graded Index (62.5/125) 800 - 1350 nm 0.275 +/- 0.015 62.5 +/- 2.5 microns 125 +/- 1 microns Bandwidth @ 850nm Greater than 200 MHz-Km Bandwidth @ 1300 nm Greater than 800 MHz-Km Better than 3.2 db per Km Better than 1.0 1.496 @ 850nm 1.491 @ 1300nm

#### 6. Contents of Tester MCS04

**6.1 Laser Diode Unit:** This comprises a laser diode operating at 650nm with a peak output power of 3 mw. Optical power from the laser is coupled to the fiber through the ceramic 2.5 mm ferrule. The ferrule is the precision cylindrical ceramic part at the tip of the cable. Optical power coupled into the fiber may be adjusted by positioning the ferrule suitably, back or forth. Please note that cable connector fits into the notch provided on the Laser Diode Unit, when the ferrule is fully inserted. Rotating the ferrule will result in variation of the pattern observed on the display screen.

The output power of the laser is settable through the potentiometer provided. One needs to set the laser output to the desired intensity. As the coupled power into the 9/125 single mode fiber will be very small, the lighting in the room may have to be dimmed, to observe the patterns with clarity.


An external 9 Vdc/500 mA power supply is required.



**Tester MCS04: Study of Mode Characteristics of Glass Optical Fibers** 

**6.2 Detector Screen:** This comprises a fixed a white screen, 54 mm wide and a cable holder on a solid metal frame. The position of the cable holder is settable at fixed distances of 30, 45, 60, 75 and 90 mm from the screen. Two M3 screws facilitate mounting of the cable holder on the frame. The actual distance between the cable end and screen needs to be measured accurately.

The cable holder has an ST adaptor that provides for holding the cable in position. Provision is also made for an



SMA cable adaptor, in case one desires to study PMMA cables with SMA connectors too.

Diameter of light spots on the scale may be measured directly employing a suitable scale. Alternately, the spot diameter may be also be computed by

capturing the image on a digital camera. The width of the screen (54 mm) or any other marking on the screen, also captured on the digital image, will serve as the reference distance in computing the actual spot diameter




- 6.3 One ST-ST 9/125 single mode step-index glass fiber patch cord.
- 6.4 One ST-ST 62.5/125 multimode graded-index glass fiber patch cord.
- 6.5 One ST-ST 50/125 multimode graded-index glass fiber patch cord (optional)


#### 7. Full Set-up to Study Single Mode Fiber Characteristics

The complete experimental set-up to study the mode characteristics of a single mode glass fiber is shown below and is self explanatory. The single mode fiber is connected to the laser diode unit and the detector screen. Laser intensity is adjusted by setting the potentiometer. Power coupled into the fiber is set by adjusting the depth of the ferrule in the laser unit. Slight tilt of the ferule may also be required for some cables

The pattern on the screen may be varied by gently distorting the cable in the midsection, to change path length of the light propagating through the core. This has to be done carefully, without disturbing the power coupling. The spot on the screen will be more distinct if the lights are dimmed in the laboratory.



#### 8. Close-up of Laser Diode Unit




9. Close-up of Detector Screen




Tester MCS04: Study of Mode Characteristics of Glass Optical Fibers


### 10. LP01 Intensity Pattern Observed on the Screen




11. LP11 Intensity Pattern Observed on the Screen



### 12. LP11 Intensity Pattern Observed on the Screen



### 13. LP22 Intensity Pattern Observed on the Screen



#### 14. NA Computation from Screen Shot of LP01 Pattern

Numerical Aperture of the 9/125 single mode step index fiber is computed using the screen shot captured during the study of the pattern. The method is described below:

Step1 Open the screen image using a suitable software (we used Inkscape) to measure the spot diameter in relation to the screen width. The actual screen width is 54mm. In the case shown below, the spot diameter is (34X54)/166.5= 11.0 mm The spot radius is 5.5 mm.

Step2: The spacing between the fiber end and the screen is fixed at 50mm.

Acceptance angle  $\theta_c$  is given by 2 tan<sup>-1</sup>(5.5/50) = 12.6 degrees

Step3: NA is given by Sin  $\theta_c/2 = 0.11$ . This matches the manufacturer's data.

Step4: V, N and other parameters of the optical fiber cable may be calculated from using the manufacturer's data and the NA obtained above. It may be noted that V is greater than 2.405 and hence at 650nm this fiber is a multimode fiber.


NA MEASUREME	NT : SI SINGLEMODE	NOTE:
34mm		
166.5mm		
MEASURED NA=0.102		

Note: This example refers to a screen-fiber end spacing of 50 mm


#### 15. Full Set-up to Study Multimode GI Fiber Characteristics

The complete experimental set-up to study the mode characteristics of a multimode graded index glass fiber is shown below and is self explanatory. The multimode graded index fiber is connected to the laser diode unit and the detector screen. Laser intensity is adjusted by setting the potentiometer. Power coupled into the fiber is set by adjusting the depth of the ferrule in the laser unit. Slight tilt of the ferule may also be required for some cables

The pattern on the screen may be varied by gently distorting the cable in the midsection, to change path length of the light propagating through the core. This has to be done carefully, without disturbing the power coupling. The spot on the screen will be more distinct if the lights are dimmed in the laboratory.




#### 16. Close-up of Laser Diode Unit




#### 17. Close-up of Detector Screen Diode Unit



#### 18. Intensity Patterns Observed on the Screen





#### **19. NA Computation from Screen Shot of Multimode Fiber Spot**


Numerical Aperture of the 62.5 /125 multimode graded index fiber is computed using the screen shot captured during the study of the pattern. The method is described below:

Step1 Open the screen image using a suitable software (we used Inkscape) to measure the spot diameter in relation to the screen width. The actual screen width is 54mm. In the case shown below, the spot diameter is (116X54)/219= 28.6 mm The spot radius is 14.3 mm

Step2: The spacing between the fiber end and the screen is fixed at 50mm. Acceptance angle  $\theta_c$  is given by 2 tan<sup>-1</sup>(14.3/50) = 32 degrees

Step3: NA is given by Sin  $\theta_c/2 = 0.275$ . This matches the manufacturer's data.

Step4: V, N and other parameters of the optical fiber cable may be calculated from using the manufacturer's data and the NA obtained above.



Note: This example refers to a measurement done with 50 mm spacing between screen and the fiber tip